

Nucleoside kinases, rate-limiting step of nucleoside analogues activation

Nucleoside analogues have proven to be a highly successful class of anti-cancer and anti-viral drugs. The therapeutic efficacy of nucleoside analogues is dependent of their intracellular phosphorylation. Two cellular nucleoside kinases, deoxycytidine kinase (dCK) and UMP-CMP kinase (CMK) are critical for phosphorylation of cytidine analogues. These kinases provide two first steps of activation of highly effective anti-cancer and anti-viral drugs, such as 1- β -D-arabinofuranosylcytosine (araC, **aracytidine**), 2',2'difluorodeoxycytidine (dFdC, **gemcitabine**), β -D-2'3'-dideoxycytidine (ddC). Both kinases phosphorylate unnatural L-nucleosides (e.g., β -L-2'3'-dideoxy-3'thiacytidine, L-SSdC, 3-TC or **lamividune**). Kinetic constants of araC, dFdC and 3TC phosphorylation by recombinant dCK and UMP-CMPK have been published. The comparison of phosphorylation properties of new nucleoside analogues with those of known drugs provides the rational basis for selection of analogues of better therapeutic potential.

To characterize the phosphorylation properties of new nucleoside analogues, **NovoCIB** has developed human recombinant dCK and human recombinant CMK nucleoside phosphorylation assays. As shown in Table 1, CMK assay must be performed with monophosphate forms of nucleoside analogues and requires preliminary phosphorylation of nucleoside analogues and their purification. To circumvent this time-consuming step, **NovoCIB** has developed a coupled dCK-CMK nucleoside phosphorylation assay that delivers in one step the critical information on both dCK and CMK substrate properties of nucleoside analogue.

Ribavirin (1-β-D-ribofuranosyl-1,2,4-triazole-3-carboxamide) is a purine nucleoside analogue with a broad-spectrum antiviral activity. Since the 1970's¹, it is known that the initial step of ribavirin phosphorylation is provided by adenosine kinase. Recently it has been demonstrated that cytosolic 5'-nucleotidase II can also phosphorylate ribavirin, that could contribute to the development of ribavirin-induced haemolytic anemia *in vivo*². **NovoCIB** has developed both human recombinant adenosine kinase and cytosolic nucleotidase II nucleoside phosphorylation assays to evaluate the properties of new ribonucleoside analogues in comparison with those of ribavirin.

	dCK assay	CMK assay	Coupled dCK-CMK assay	AK assay	5'cN-II assay
Natural substrates	Deoxyadenosine Deoxyguanosine Deoxycytidine Cytidine	dCMP, UMP CMP	Deoxycytidine Cytidine	Adenosine Inosine	Deoxyinosine Inosine
Nucleoside analogues substrates	Cladribine, fludarabine Gemcitabine, Lamivudine, Aracytidine Fluorodeoxyuridine	Monophosphate forms of cytidine analogues dFdCMP, araCMP, 3TCMP Adefovir (9-(2- phosphonomethoxyethyl) adenine	Gemcitabine Aracytidine Lamivudine Fluodeoxyuridine	Ribavirin Tubercidin Mizoribin	Dideoxyinosine Ribavirin Acyclovir

Table 1. Available nucleoside kinase assays and reference nucleoside substrates

_

¹ R. C. Willis, D. A. Carson, and J. E. Seegmiller (1978) Adenosine Kinase Initiates the Major Route of Ribavirin Activation in a Cultured Human Cell Line. PNAS USA, 75: 3042-3044

² Wu JZ, Larson G, Walker H, Shim JH, Hong Z. Phosphorylation of ribavirin and viramidine by adenosine kinase and cytosolic 5'-nucleotidase II: Implications for ribavirin metabolism in erythrocytes (2005) Antimicrob Agents Chemother. 49(6):2164-71.



Ref: # E-Nov 3

Human deoxycytidine kinase (dCK) Human, recombinant expressed in E.coli E.C. 2.7.1.74

Description

NOVOCIB's human deoxycytidine kinase (dCK) is a recombinant protein of ca.33kDa cloned by RT-PCR amplification of mRNA extracted from human hepatoma cells and expressed in E.coli.


Human deoxycytidine kinase plays a key role in the salvage pathway of deoxynucleotides synthesis providing resting cells with deoxynucleotides for DNA repair and mitochondrial DNA synthesis. The enzyme has a broad substrate specificity and provides the phosphorylation of both purine and pyrimidine deoxynucleosides (e.g. deoxyadenosine (dA), deoxyguanosine (dG)) and deoxycytidine (dC) and pyrimidine ribonucleoside, cytidine (C)). The enzyme can utilize both ATP and UTP as phosphate donor with UTP as a preferred substrate.


Deoxycytidine kinase is responsible for the phosphorylation and activation of numerous nucleoside analogs used to treat cancer (e.g. cytarabine, gemcitabine, cladribine and fludarabine) including nucleoside analogs of non-physiological L-chirality (e.g. 3TC, lamivudine, anti-HIV and anti-hepatitis B agent). Three-dimensional structures of dCK in complex with various pyrimidine¹ and purine^{2,3} D- and L-nucleosides have been solved providing structural basis for activation of L- and D-nucleoside analogs.

Storage: -20 ℃ in a solution containing 50 mM Tris-HCl, pH 7.6, 1 mM β-mercaptoethanol, 50% glycerol.


Unit Definition: One unit of deoxycytidine kinase converts 1.0 µmole of deoxycytidine and ATP to dCMP and ADP per minute at pH 7.6 at 37℃, as measured by a coupled PK/LDH enzyme system.

Specific Activity: 0.025 unit/mg protein. Purity: controlled by 12%AA





activity of enzvmatic human recombinant dCK was confirmed by ionpair HPLC analysis (Agilent 1100 series, Zorbax C18plus) as shown by formation dCMP and ADP (red) from deoxycytidine and ATP (blue).



Assay condition: Enzymatic activity of dCK is measured by spectrophotometric assays in a coupled lactate dehydrogenase/pyruvate kinase system. Assays were carried out at 37°C, at 50mM Tris-HCl pH7,6; 50mM KCI, 10mM MgCl2, 5mM ATP, 0,1mM NADH, 1mM phosphoenolpyruvate, 1mM DTT, PK 10U/ml, LDH 15U/ml, 0,9µM dCK. Reaction was followed in an iEMS Reader MF (Labsystems) microtiter plate reader at 340nm. Nucleosides, nucleotides, LDH and PK were purchased from Sigma-Aldrich.

Related products:

NOVOCIB has cloned and purified a panel of human recombinant nucleoside kinases and has developed a range of PRECICE® services to evaluate substrate properties of new nucleoside analogues for key cellular kinases.

- dCK nucleoside phosphorylation assay
- Coupled dCK-CMK nucleoside phosphorylation assays
- Coupled Nucleoside Kinase IMPDH II
- UMP-CMP kinase (CMK)
- Adenosine kinase (AK)
- Cytosolic 5' nucleotidase II (cN-II)
- CMK nucleotide monophosphate phosphorylation assay
- Adenosine kinase phosphorylation assay
- cN-II phosphorylation assay

NovoCIB SAS, 115 avenue Lacassagne, 69003 Lyon, France contact@novocib.com Tel / Fax +33 (0)478536395

Sabini E, Ort S, Monnerjahn C, Konrad M, Lavie A. Structure of human dCK suggests strategies to improve anticancer and antiviral therapy. (2003) Nat Struct

Sabini E, Hazra S, Ort S, Konrad M, Lavie A. Structural basis for substrate promiscuity of dCK. (2008) J Mol Biol. 378(3):607-21

³ Sabini E, Hazra S, Konrad M, Burley SK, Lavie A. (2007) Structural basis for activation of the therapeutic L-nucleoside analogs 3TC and troxacitabine by human deoxycytidine kinase. *Nucleic Acids Res.* 35(1):186-92



Ref: # IVS-Nov 3

dCK nucleoside phosphorylation assay

IMPORTANT: Client-specified alterations can be accommodated


Aim: Characterization of substrate properties (Km and Vmax) of new nucleoside analogues for human deoxycytidine kinase in comparison with the properties of known nucleoside analogues (e.g. aracytidine, gemcitabine, cladribine and lamivudine).

	Novocib*			Published data			
Substrate	Km,µM	Vmax, µmol/mg/min	Relative Vmax, % of dCR	Km,µM	Vmax, µmol/mg/mi n	Ref.	
	0,577	0,026	100	0,16	0,033	Recombinant Johansson Karlsson 1995 ¹	
Deoxycytidine				1,3	0,069	Recombinan Usova & Eriksson, 1997 t ²	
				0,57	0,004	Partially purified Someya H et al 2003 ³	
Gemcitabine	42,71	0,325	1250				
Deoxyadenosine	150,5	1,08	4153	115		Recombinant Sabini E et al 2008 ⁴	
Deoxyaderiosirie	130,3	1,00	4133	480	1,5	Recombinant Johansson Karlsson 1995	
Aracytidine	6,81	0,224	862	15	0,009	Partially purified Someya H et al 2003	
Cladribine	56,5 0,285	0.285	1096	89	0,126	Recombinant Usova & Eriksson, 1997	
		3,200		24	0,76	Recombinant Johansson Karlsson 1995	

Enzyme: The dCK used in the assays is a human recombinant dCK, cloned from human cells, expressed in E. coli, produced and purified by NOVOCIB (see sheet # E-Nov 3 for further information). The enzyme purity is controlled by SDS-PAGE. Protein concentration is measured by Bradford method (Bio-Rad). dCK enzymatic activity (≥ 0.025 unit/mg protein) is systematically controlled before performing any assay.

Kinetics Analysis: Enzymatic activity of deoxycytidine kinase with particular nucleoside substrate is measured continuously by spectrophotometric assays in a coupled lactate dehydrogenase/pyruvate kinase system. Assays are carried out at 37°C, at 50mM Tris-HCl pH7,6; 50 mM KCI, 10mM MgCI2, 5mM ATP, 0,1mM NADH, 1mM phosphoenolpyruvate, 1mM DTT, PK 10U/ml, LDH 15U/ml, 0,9µM dCK. The nucleosides, nucleotides, LDH and PK are purchased from Sigma-Aldrich. Reaction is followed in an iEMS Reader MF (Labsystems) microtiter plate reader at 340nm. Assays are performed in duplicate (2 wells per compound and per concentration). Triplicates are available upon request. Km and Vmax are calculated from spectroscopic data using Michaelis-Menten equation.

A confirmation by HPLC analysis of formation of monophosphorylated forms is available upon request.



Related products:

NOVOCIB has cloned and purified a panel of human recombinant nucleoside kinases and has developed a range of PRECICE® services to evaluate substrate properties of new nucleoside analogues for key cellular kinases.

- Deoxycytidine kinase (dCK)
- Coupled dCK-CMK nucleoside phosphorylation assays
- Coupled Nucleoside Kinase IMPDH II
- UMP-CMP kinase (CMK)
- Adenosine kinase (AK)
- Cytosolic 5' nucleotidase II (cN-II)
- CMK nucleotide monophosphate phosphorylation assay
- Adenosine kinase phosphorylation assay
- cN-II phosphorylation assay

¹ M. Johansson and A. Karlsson (1995): Differences in kinetic properties of pure recombinant human and mouse deoxycytidine kinase Biochem. Pharmacol. 50(2), 163-

^{168 &}lt;sup>2</sup> E. V. Usova and S. Eriksson (1997) The effects of high salt concentrations on the regulation of the substrate specificity of human recombinant deoxycytidine kinase Eur. J. Biochem. 248(3), 762-766

H. Someya et al. (2003) Phosphorylation of 4'-thio-beta-D-arabinofuranosylcytosine and its analogs by human deoxycytidine kinase J. Pharmacol. Exp. Ther. 304(3),

<sup>1314-1322

&</sup>lt;sup>4</sup> E. Sabini *et al.* (2008) **Structural basis for substrate promiscuity of dCK** *J. Mol. Biol.* 378(3), 607-621



PRECICE® Services Information sheet Ref: # E-Nov 4

UMP-CMP kinase (CMK) Human, recombinant expressed in E.coli

E.C. 2.7.4.14

Synonyms: cytidylate kinase, deoxycytidylate kinase, deoxycytidine monophosphokinase, dCMP kinase, cytidine monophosphate kinase, CMP kinase (CMK, CMPK), uridine monophosphate kinase (UMK, UMPK), uridine monophosphate/cytidine monophosphate kinase, UMP/CMP kinase (UMP/CMPK), CTP:CMP phosphotranferase, ATP:UMP-CMP phosphotransferase, pyrimidine nucleoside monophosphate kinase (YMPK)

Description

NOVOCIB's Human UMP-CMP kinase (CMK) is a recombinant protein of ca. 27kDa (full length 228-aa form¹) cloned by RT-PCR amplification of mRNA extracted from Huh7 cells (human hepatoma) and expressed in E.coli.


UMP-CMP kinase plays a critical role in supplying cells with nucleotides by catalysing the phosphorylation of CMP, UMP and dCMP to their respective diphosphates. CMK plays also an important role in the activation of cytidine analogues, aracytidine and gemcitabine, a mainstay of leukaemia and lymphoma therapy². CMK has a remarkable ability of to phosphorylate L-nucleotides from their monophosphate to diphosphate forms³ as shown for β-L-2'3'-dideoxy-3'thiacytidine (L-SSdC, 3-TC or lamividune), an anti-HIV and anti-hepatitis B drug.


Crystal structure of open form of human UMP-CMP kinase has been solved recently⁴. These data, together with the homology model of enzyme in closed state, provides structural basis for understanding the substrate specificity of the enzyme and helps to design new nucleoside analogues of higher phosphorylation efficiency.

Storage: -20 ℃ in a solution containing 150mM KCI, 50mM Tris-Hcl, pH7,5, 2mM βmercaptotethanol, 50% glycerol.


Unit Definition: One unit of UMP-CMP kinase converts 1.0 µmole of UMP and ATP to UDP and ADP per minute at pH 7.6 at 25℃, using a coupled enzyme system with PK/LDH.

Specific activity: ≥0,150U/mg Purity: controlled by SDS-PAGE





UMP-CMP kinase was confirmed by HPLC analysis as illustrated by CDP and ADP formation (red) from CMP and ATP (blue).



Assay condition: Enzymatic activity of UMP-CMP kinase is measured by continuous spectrophotometric assays in a coupled lactate dehydrogenase/pyruvate kinase system. Assays are carried out at 37°C, at 50mM Tris-HCl pH7,6; 50mM KCl, 10mM MgCl₂, 5mM ATP, 0,1mM NADH, 1mM phosphoenolpyruvate, 1mM DTT, PK 10U/ml, LDH 15U/ml, 380nM CMK. Reaction is followed in an iEMS Reader MF (Labsystems, Finland) microtiter plate reader at 340nm.

Related products:

NOVOCIB has cloned and purified a panel of human recombinant nucleoside kinases and has developed a range of PRECICE® services to evaluate substrate properties of new nucleoside analogues for key cellular kinases.

- UMP-CMP kinase (CMK) nucleoside phosphorylation assay
- Coupled dCK-CMK nucleoside phosphorylation assays
- Deoxycytidine kinase (dCK)
- Adenosine kinase (AK)
- Cytosolic 5' nucleotidase II (cN-II)
- dCK nucleoside phosphorylation assay
- Adenosine kinase phosphorylation assay
- cN-II phosphorylation assay
- Coupled Nucleoside Kinase IMPDH II

NovoCIB SAS, 115 avenue Lacassagne, 69003 Lyon, France Tel / Fax +33 (0)478536395

Jieh-Yuan Liou, Ginger E. Dutschman, Wing Lam, Zaoli Jiang and Yung-Chi Cheng (March 2002) Characterization of Human UMP/CMP Kinase and Its Phosphorylation of D- and L-Form Deoxycytidine Analogue Monophosphates Cancer Research 62, 1624-1631
 Van Rompay AR, Johansson M, and Karlsson A (September 1999) Phosphorylation of Deoxycytidine Analog Monophosphates by UMP-CMP Kinase: Molecular

Characterization of the Human Enzyme Mol Pharmacol 56 (3), 562-569

Claudia Pasti, Sarah Gallois-Montbrun, Hélène Munier-Lehmann, Michel Veron, Anne-Marie Gilles and Dominique Deville-Bonne (Mar2003) Reaction of human UMP-CMP

kinase with natural and analog substrates. European Journal of Biochemistry, 270 (8), 1784 - 1790

Segura-Peña D, Sekulic N, Ort S, Konrad M, Lavie A. (2004) Substrate-induced conformational changes in human UMP/CMP kinase. J Biol Chem. 279(32):33882-33889



CMK nucleotide monophosphate phosphorylation assay

IMPORTANT: Client-specified alterations can be accommodated.


Aim: Characterization of substrate properties (Km and Vmax) of monophosphate forms of new nucleoside analogues for human CMK in comparison with monophosphate forms of natural nucleosides or reference nucleoside analogues.

Enzyme: The enzyme used in the assays is a human recombinant CMK, cloned from human cells, expressed in *E. coli*, produced and purified by **NOVOCIB** (see sheet # E-Nov 4 for further information). The enzyme purity is controlled by SDS-PAGE. Protein concentration is measured by Bradford method (Bio-Rad). CMK enzymatic activity (≥ 0.150 unit/mg protein) is systematically controlled before performing any assay.

			NovoCIB*	Published data ⁶	
	Km, μM	Vmax nmol/mg/min	Relative Vmax, % of CMP	Km, μM	Vmax, nmol/mg/min
CMP	17.9	130.07	100	20	350
UMP	392	307.16	236	45	350
dCMP	1334	297.65	228	900	200

Kinetics Analysis: Substrate properties of a particular nucleoside monophosphate for CMK are evaluated in a continuous LDH/PK spectrophotometric assay. The assays are carried out at 37°C, at 50mM Tris-HCl pH 7.6; 50mM KCl, 10mM MgCl₂, 5mM ATP, 0.1mM NADH, 1mM phosphoenolpyruvate, 1mM DTT, PK 10U/ml, LDH 15U/ml, 380nM CMK. Nucleosides, nucleotides, LDH and PK are purchased from Sigma-Aldrich. Reaction is followed in an iEMS Reader MF (Labsystems) microtiter plate reader at 340nm. Assays are performed in duplicate (2 wells per compound and per concentration). Triplicates are available upon request. Km and Vmax are calculated from spectroscopic data using Michaelis-Menten equation.

A confirmation by HPLC analysis of the formation of monophosphorylated forms is available upon request.



Related products:

NOVOCIB has cloned and purified a panel of human recombinant nucleoside kinases and has developed a range of PRECICE® services to evaluate substrate properties of new nucleoside analogues for key cellular kinases.

- UMP-CMP kinase (CMK)
- Coupled dCK-CMK nucleoside phosphorylation assays
- Deoxycytidine kinase (dCK)
- Adenosine kinase (AK)
- Cytosolic 5' nucleotidase II (cN-II)
- · dCK nucleoside phosphorylation assay
- Adenosine kinase phosphorylation assay
- cN-II phosphorylation assay
- Coupled Nucleoside Kinase IMPDH II


References

Chih-Hung Hsu, Jieh-Yuan Liou, Ginger E. Dutschman, and Yung-Chi Cheng (2005) Phosphorylation of Cytidine, Deoxycytidine, and Their Analog Monophosphates by Human UMP/CMP Kinase is Differentially Regulated by ATP and Magnesium *Mol Pharmacol* 67:806-814
Topalis D, Kumamoto H, Amyav eValasco MF, Dugué L, Haouz A, Alexandre JA, Gallois-Montbrun S, Alzari PM, Pochet S, Agrofoglio LA, Deville-Bonne D. (Jul 2007) Nucleotide binding to human UMP-CMP kinase using fluorescent derivatives – a screening based on affinity for the UMP-CMP binding site. *FEBS J.* 274(14):3704-14





Coupled dCK-CMK nucleoside phosphorylation assays

Aim: Coupled dCK-CMK nucleoside phosphorylation assay is a cost-effective rapid assay that delivers in one step the critical information on both dCK and CMK substrate properties of nucleoside analogue.



Enzymes: The enzymes used in this assay are human recombinant dCK and human recombinant CMK, cloned by NovoCIB from human cells, expressed in *E. coli*, and produced and purified by **NOVOCIB** (see sheet # E-Nov 3 and # E-Nov 4 for detailed information).





Assay condition:

Substrate properties of nucleoside analogue for dCK and CMK kinases evaluated in a two-step spectrophotometric assay carried out at 37°C, at 50mM Tris-HCl pH7,6; 50mM KCl, 10mM MgCl $_2$, 5mM ATP, 0.1mM NADH, phosphoenolpyruvate, 1mM DTT, PK 10U/ml, LDH 15U/ml. The initial phosphorylation is started by addition of dCK (1µM) and the reaction is followed spectropho-metrically 340nm during 90 min followed by addition of CMK (0.3µM). The changes in absorbance at 340nm are used to calculate both the initial rate of reactions and the concentration of nucleoside mono-phosphate formed.

Method validation:

The phosphorylation kinetic of CMP by recombinant CMK have been measured in two independent approaches. In first one, CMK Km for CMP was studied directly with CMP substrate (grey), and in second one CMK Km for CMP was measured indirectly in coupled dCK-CMK assay (red) usine cytidine as a substrate. As shown on left, coupled dCK-CMK assay produces results which are highly similar to those obtained from a direct CMK assay.

Related products:

NOVOCIB has cloned and purified a panel of human recombinant nucleoside kinases and has developed a range of PRECICE® services to evaluate substrate properties of new nucleoside analogues for key cellular kinases.

- UMP-CMP kinase (CMK) nucleoside phosphorylation assay
- dCK nucleoside phosphorylation assay
- UMP-CMP kinase (CMK)
- Deoxycytidine kinase (dCK)
- Adenosine kinase (AK)
- Cytosolic 5' nucleotidase II (cN-II)
- Adenosine kinase phosphorylation assay
- cN-II phosphorylation assay
- Coupled Nucleoside Kinase IMPDH II



Ref: # E-Nov 5

Human adenosine kinase (AK) Human, recombinant expressed in E.coli EC 2.7.1.20

Synonyms: ADK, Adenosine 5'-phosphotransferase

Description

NOVOCIB's human adenosine kinase (AK) is a recombinant protein of ca.39kDa (345-aa short form^{1, 2}) cloned by RT-PCR amplification of mRNA extracted from human hepatoma cells and expressed in E.coli. The sequence of the cloned AK (GenBank accession number U50196) was confirmed by DNA sequencing (100% identity).

Adenosine kinase is a ubiquitous enzyme that catalyzes the transfer of γ-phosphate from ATP to 5' hydroxyl of adenosine generating AMP and ADP. Adenosine (AR) is an important modulator of central nervous system functions with a half-life of seconds. Facilitated diffusion of adenosine across the cell membrane closely couples adenosine concentrations in the intracellular and extracellular compartments. Inhibition of adenosine kinase results in selective increase of local adenosine concentrations and reduced seizure susceptibility and nociception in vivo3. Adenosine kinase is an attractive and experimentally validated target for the development of new analgesic and anti-inflammatory agents4. In addition, AK recently has emerged as a novel target to predict and to prevent epileptogenesis^{5, 6}. The X-ray crystallographic structure of human AK has been described⁷ and provides structural basis for rational design and optimisation of new AK inhibitors.


In addition, this enzyme is responsible for the phosphorylation and consequent clinical activity of several therapeutically useful nucleosides, including the antiviral drug ribavirin⁸, immunosuppressive drug mizoribine⁹ and anticancer C-nucleoside, tiazofurin¹

Storage: -20 ℃ in a solution containing 50 mM Tris-HCl, pH 7.6, 1 mM β-mercaptoethanol, 50% glycerol.

Unit Definition: One unit of adenosine kinase converts 1.0 µmole of adenosine and ATP to AMP and ADP per minute

at pH 7.6 at 30℃, as measured by a coupled PK/LDH enzyme system.

Specific Activity: \geq 0.030 unit/mg protein. Purity: controlled by 10% AA SDS-PAGE.



Assay condition: Enzymatic activity of adenosine kinase with particular nucleoside substrate is measured by spectrophotometric assays in a Deoxyadenosine Km = 295μ coupled lactate dehydrogenase / pyruvate kinase system. Assays were carried out at 37℃, at 50mM Tris-HCI pH7,6; 50mM KCI, 5mM MgCI2, 2,5mM ATP, 0,1mM NADH, 1mM phosphoenolpyruvate, 1mM DTT, PK 5U/ml, LDH 5U/ml. Reaction was followed in an iEMS Reader MF (Labsystems) microtiter plate reader at 340nm. Nucleosides, nucleotides, 12LDH and PK were purchased from Sigma-Aldrich.

Related products:

NOVOCIB has cloned and purified a panel of human recombinant nucleoside kinases and has developed a range of PRECICE® services to evaluate substrate properties of new nucleoside analogues for key cellular kinases.

- Adenosine kinase phosphorylation assay
- Coupled Nucleoside Kinase IMPDH II
- Deoxycytidine kinase (dCK)
- UMP-CMP kinase (CMK)
- Cytosolic 5' nucleotidase II (cN-II)
- CMK nucleotide monophosphate phosphorylation assay
- dCK nucleoside phosphorylation assay
- Coupled dCK-CMK nucleoside phosphorylation assays
- cN-II phosphorylation assay

¹⁰ Saunders PP, Spindler CD, Tan MT, Alvarez E, Robins RK **Tiazofurin Is Phosphorylated by Three Enzymes from Chinese Hamster Ovary Cells** (1990) Cancer Research 50, 5269-5274

NovoCIB SAS, 115 avenue Lacassagne, 69003 Lyon, France contact@novocib.com Tel / Fax +33 (0)478536395

^{1.} Spychala J, Datta NS, Takabayashi K, Datta M, Fox IH, Gribbins T, Mitchell BS Cloning of human adenosine kinase cDNA: Sequence similarity to microbial ribokinases and fructokinases (1996) *Proc. Natl. Acad. Sci. USA* 93, pp. 1232-1237
² Sahin B, Kansy JW, Nairn AC, Spychala J, Ealick SE, Fienberg AA, Greene RW and Bibb JA. **Molecular characterization of recombinant mouse adenosine kinase and**

evaluation as a target for protein phosphorylation Eur. J. Biochem. 271, 3547–3555 (2004)

³ Jarvis MF, Yu H, Kohlhaas K, Alexander K, Lee CH, Jiang M, Bhagwat SS, Williams M, Kowaluk EA ABT-702 (4-Amino-5-(3-bromophenyl)-7-(6-morpholinopyridin-3-yl)pyrido[2,3-d]pyrimidine), a Novel Orally Effective Adenosine Kinase Inhibitor with Analgesic and Anti-Inflammatory Properties: I. In Vitro Characterization and Acute Antinociceptive Effects in the Mouse Journal of Pharmacology and Experimental Therapeutics 295:1156–1164, 2000

McGaraughty S, Cowart M, Jarvis MF, Berman RF, Anticonvulsant and antinociceptive actions of novel adenosine kinase inhibitors Curr Top Med Chem. 2005;5(1):43-

Fedele DE, Gouder N, Guttinger M, Gabernet L, Scheurer L, Rulicke T, Crestani F, Boison D **Astrogliosis in epilepsy leads to overexpression of adenosine kinase, resulting in seizure aggravation** (2005) *Brain*, 128, 2383–2395

Detlev Boison **The adenosine kinase hypothesis of epileptogenesis** *Progress in Neurobiology* 84 (2008) 249–262

Mathews, I.I., Erion, M.D. & Ealick, S.E. Structure of human adenosine kinase at 1.5 A resolution. (1998) Biochemistry 37, 15607–15620 Willis RC, Carson DA, Seegmiller JE. Adenosine Kinase Initiates the Major Route of Ribavirin Activation in a Cultured Human Cell (1978) PNAS 1978;75;3042-3044

⁹ Miller RL, Adamczyk DL, Miller WH, Koszalka GW, Rideout JL, Beacham LM, Chao EY, Haggerty JJ, Krenitsky TA, Elion, GB Adenosine kinase from rabbit liver II Substrate and inhibitor specificity. (1979) *J. Biol. Chem.* 254, 2346-2352






AK nucleoside phosphorylation assay

IMPORTANT: Client-specified alterations can be accommodated.

Aim: Characterization of substrate properties (Km and Vmax) of new nucleoside analogues for human adenosine kinase in comparison with properties of known nucleoside analogues (e.g. ribavirine, tubercidine or mizoribine).

	Novocib		Published		
Substrate	Km (µM)	Kcat (min ⁻¹)	Km (µM)	Kcat (min ⁻¹)	Ref
Adenosine	11	1.5	3.2	13	1
			0.150		2
Ribavirine	328	1,9	540	1,8	1
Deoxyadenosine	295	3,4	360		2
Tubercidine	12	2,2			
Inosine	1758	2,6			


Enzyme: The AK used in the assays is a human recombinant AK, cloned from human cells, expressed in E. coli, produced and purified by NOVOCIB (see sheet # E-Nov5 for further information). The enzyme purity is controlled by SDS-PAGE. Protein concentration is measured by Bradford method (Bio-Rad). AK enzymatic activity (≥ 0.030 unit/mg protein) is systematically controlled before performing any assay.



The phosphorylation of ribavirin bv adenosine kinase was confirmed by as analysis illustrated Ribavirine-MP formation (red) ribavirine (blue).

Related products:

NOVOCIB has cloned and purified a panel of human recombinant nucleoside kinases and has developed a range of PRECICE® services to evaluate substrate properties of new nucleoside analogues for key cellular kinases.



Kinetics Analysis: Enzymatic activity of adenosine kinase with particular nucleoside substrate is measured continuously by spectrophotometric assays in a coupled lactate dehydrogenase/pyruvate kinase system. Assays are carried out at 37°C, at 50mM Tris-HCl pH7,6; 50mM KCl, 5mM MgCl2, 2,5mM ATP, 0,1mM NADH, 1mM phosphoenolpyruvate, 1mM DTT, PK-LDH (5U/ml each), 0,85 μ M AK. The nucleosides, nucleotides, LDH and PK are purchased from Sigma-Aldrich. Reaction is followed in an iEMS Reader MF (Labsystems) microtiter plate reader at 340nm. Assays are performed in duplicate (2 wells per compound and per concentration). Triplicates are available upon request. Km and Vmax are calculated from spectroscopic data using Michaelis-Menten equation.

A confirmation by HPLC analysis of formation of monophosphorylated forms is available upon request.

- Adenosine kinase
- Coupled Nucleoside Kinase IMPDH II
- Deoxycytidine kinase (dCK)
- UMP-CMP kinase (CMK)
- Cytosolic 5' nucleotidase II (cN-II)
- CMK nucleotide monophosphate phosphorylation assay
- dCK nucleoside phosphorylation assay
- Coupled dCK-CMK nucleoside phosphorylation assays
- cN-II phosphorylation assay

^{1.} Wu JZ, Larson G, Walker H, Shim JH, Hong Z. Phosphorylation of ribavirin and viramidine by adenosine kinase and cytosolic 5'-nucleotidase II: Implications for ribavirin metabolism in erythrocytes. (2005) Antimicrob Agents Chemother. 49(6):2164-71
2. Yamada, Y.; Goto, H.; Ogasawara, N. Adenosine kinase from human liver (1981) Biochim. Biophys. Act, 660, 36-43



Ref: E-Nov 6

Human cytosolic 5'-nucleotidase II (cN-II) Human, recombinant expressed in E.coli

EC 3.1.3.5

Synonyms: 5'-nucleotidase/phosphotransferase, cytosolic High Km 5'-nucleotidase (hkm-NT), cytosolic purine 5'-nucleotidase (purine 5'-NT), IMP/GMPspecific 5'-nucleotidase (IMP/GMP-specific 5'-NT)

Description

NOVOCIB's human cytosolic IMP/GMP specific 5'-nucleotidase/phosphotransferase II (cN-II) is a recombinant protein of ca. 65kDa cloned by RT-PCR amplification of mRNA extracted from human hepatoma cells and expressed in E.coli. The sequence of the cloned NT5C2 gene (GenBank accession number P49902) was confirmed by DNA sequencing (100% identity).

Cytosolic 5'-nucleotidase II is one of the seven known mammalian nucleotidases¹ that specifically catalyzes the dephosphorylation of 6-hydroxypurine nucleoside 5'-monophosphates (IMP, dIMP, dGMP) and regulates cellular pool of IMP and GMP^{2, 3}. The enzyme also acts as a phosphotransferase catalyzing the transfer of a phosphate from nucleoside monophosphate to a nucleoside acceptor - preferentially inosine and deoxyinosine. Unlike the other 5'-nucleotidases, cN-II is allosterically regulated by adenine/guanine nucleotides and 2,3-biphosphoglycerate⁴.

In addition, cytosolic 5'-nucleotidase II phosphorylates anti-viral and anti-tumour nucleoside analogues such as 2'3'-dideoxyinosine, carbovir⁵, acyclovir⁶ and ribavirin⁷.

Storage: -20 ℃ in a solution containing 50 mM Tris-HCl, pH 7.6, 2 mM β-mercaptoethanol, 50% glycerol.



Unit Definition: One unit of 5'nucleotidase converts 1.0 µmole of IMP to inosine per minute at pH 7.6 at 37℃, as measured by a coupled PNP/XO enzyme system in the presence of 20mM MgCl2, 5mM DTT, 500µM KH₂PO₄, and 1,25mM IMP.

Specific Activity:

≥ 0.150 unit/mg protein.

Purity: controlled by 10% AA SDS-

PAGÉ.

Vi, µM/min 2.5 IMP Km = 0,945mM Vmax = 0,465U/mg with 1.5mMdATP IMP Km = 0.34mM Vmax = 0.518U/mg IMP, mM

5'-nucleotidase assay condition: 5'-nucleotidase activity of cN-II is followed in an irreversible spectrophotometric assay using coupled purine nucleoside phosphorylase - xanthine oxidase system (2,5mU/ml each). Assays were carried out at 37°C, at 50mM Tris-HCl pH7,6; 100mM KCl, 20mM Mg Cl₂, 500µM KH₂PO₄, 5mM DTT, 119nM cN-II and various concentration of IMP. Reaction is followed at 295nm. The IMP is purchased from MP Biochemicals and XO from Sigma-Aldrich. PNP is produced and purified by NovoCIB (ref # E-Nov2).

Related products:

NOVOCIB has cloned and purified a panel of human recombinant nucleoside kinases and has developed a range of PRECICE® substrate services to evaluate new nucleoside analogues for key cellular kinases.

- cN-II phosphorylation assay
- Coupled Nucleoside Kinase IMPDH II
- · Adenosine kinase
- Deoxycytidine kinase (dCK)
- UMP-CMP kinase (CMK)
- dCK nucleoside phosphorylation assay
- CMK nucleotide monophosphate phosphorylation assay
- Coupled dCK-CMK nucleoside phosphorylation assays
- Adenosine kinase nucleoside phosphorylation assays

NovoCIB SAS, 115 avenue Lacassagne, 69003 Lyon, France contact@novocib.com Tel / Fax +33 (0)478536395

Bianchi V, Spychala J. Mammalian 5'-nucleotidases (2003) J. Biol Chem 278(47): 46195-46198

² Allegrini S, Pesi R, Tozzi MG, Fiol CJ, Johnson RB, Eriksson S. **Bovine cytosolic IMP/GMP-specific 5'-nucleotidase: cloning and expression of active enzyme in Escherichia coli.** (1997) *Biochem J.* 328:483-7.

lpata PL, Tozzi MG Recent advances in structure and function of cytosolic IMP-GMP specific 5'-nucleotidase II (cN-II) (2006) Purinergic Signal. 2(4):669-75 Spychala J, Madrid-Marina V, Fox IH High Km Soluble 5'-nucleotidase from human placenta (1988) J Biol Chem 263(35): 18759-18765

Johnson MA, Fridland A. Phosphorylation of 2',3'-dideoxyinosine by cytosolic 5'-nucleotidase of human lymphoid cells. (1989) Mol Pharmacol. 36(2):291-5 Keller PM, McKee SA, Fyfe JA. Cytoplasmic 5'-nucleotidase catalyzes acyclovir phosphorylation (1985) J Biol Chem. 260(15):8664-7

Wu JZ, Larson G, Walker H, Shim JH, Hong Z. Phosphorylation of ribavirin and viramidine by adenosine kinase and cytosolic 5'-nucleotidase II: Implications for ribavirin metabolism in erythrocytes (2005) Antimicrob. Agents Chemother. 49(6):2164-71



PRECICE® Services Information sheet Ref: IVS-Nov1&5

Coupled nucleoside kinase - IMPDH II assay


IMPORTANT: Client-specified alterations can be accommodated.


IMP Dehydrogenase (IMPDH, E.C. 1.1.1.205) catalyzes the pivotal step in guanine nucleotide biosynthesis. By converting inosine monophosphate (IMP) to xanthosine monophosphate (XMP), IMPDH controls the guanine nucleotide pool. A number of nucleoside analogues (e.g. ribavirin, mizoribine) are known to inhibit IMPDH after being monophosphorylated. The therapeutic consequences of IMPDH inhibition vary from different analogues - mizoribine is an immunosuppressor and ribavirin is a broad spectrum antiviral. Even if direct relationship between ribavirin antiviral action and IMPDH inhibition by ribavirin monophosphate has not been demonstrated, the depletion of cellular GTP might result in an increased frequency of ribavirin triphosphate incorporation by viral polymerase due to a lower intracellular concentration of its natural competitor.

Aim: For rapid evaluation of monophosphate forms of nucleoside analogues as IMPDH inhibitors.


	AK	5'cN-II
Natural	Adenosine	Deoxyinosine
substrates	Inosine	Inosine
Nucleoside	Ribavirin	Dideoxyinosine
analogues	Tubercidin	Ribavirin
substrates	Mizoribine	Acyclovir

Enzymes: The monophosphorylation step of nucleoside analogue is provided by one of the specific human recombinant nucleoside kinases: AK (ref. # E-Nov 5) or cN-II (ref. # E-Nov 6) produced by **NOVOCIB**. Human recombinant IMPDH 2 was cloned from human cells, expressed in *E. coli* and purified by **NOVOCIB** (see sheet # E-Nov 1 for further information). The enzyme purity is controlled by SDS-PAGE, protein concentration is measured by Bradford method (Bio-Rad). A standard operating procedure (SOP) is followed to measure enzymatic activity.





A confirmation by HPLC analysis of formation of monophosphorylated forms is available upon request.



IMPDH inhibition: Effect of monophosphorylated nucleoside analogues on human recombinant IMPDH II. Enzymatic assays performed in duplicate are carried out at 37°C in 0.1M KH $_2\text{PO}_4$ buffer pH 8.0 in the presence of 2mMDTT, 200µM NAD, 200µM IMP and 0.2 µM IMPDH II and increasing concentration of monophosphorylated nucleoside. Reaction is followed in an iEMS Reader MF (Labsystems) microtiter plate reader at 340nm.

References

1] P. Leyssen, J. Balzarini, E. De Clercq, J. Neyts (2005) The Predominant Mechanism by Which Ribavirin Exerts Its Antiviral Activity In Vitro against Flaviviruses and Paramyxoviruses Is Mediated by Inhibition of IMP Dehydrogenase J Virol 79: 1943–1947 [2] L.J. Stuyver, S. Lostia, S.E. Patterson, J.L. Clark, K. A. Watanabe, M.J. Otto and K.W. Pankiewicz (2002) Inhibitors of the IMPDH enzyme as potential antibovine viral diarrhoea virus agents Antiviral Chemistry & Chemotherapy 13:345–352

Related products:

- · dCK nucleoside phosphorylation assay
- Adenosine kinase nucleoside phosphorylation assays
- cN-II phosphorylation assay
- Deoxycytidine kinase (dCK)
- Adenosine kinase
- Cytosolic 5' nucleotidase II (cN-II)
- UMP-CMP kinase (CMK)
- · CMK nucleotide monophosphate phosphorylation assay
- Coupled dCK-CMK nucleoside phosphorylation assays