

## Human deoxycytidine kinase (dCK) Human, recombinant expressed in *E.coli* E.C. 2.7.1.74

### Description

**NOVO CIB's** human deoxycytidine kinase (dCK) is a recombinant protein of ca.33kDa cloned by RT-PCR amplification of mRNA extracted from human hepatoma cells and expressed in *E.coli*.

Human deoxycytidine kinase plays a key role in the salvage pathway of deoxynucleotides synthesis providing resting cells with deoxynucleotides for DNA repair and mitochondrial DNA synthesis. The enzyme has a broad substrate specificity and provides the phosphorylation of both purine and pyrimidine deoxynucleosides (e.g. deoxyadenosine (dA), deoxyguanosine (dG)) and deoxycytidine (dC) and pyrimidine ribonucleoside, cytidine (C)). The enzyme can utilize both ATP and UTP as phosphate donor with UTP as a preferred substrate.



Deoxycytidine kinase is responsible for the phosphorylation and activation of numerous nucleoside analogs used to treat cancer (e.g. cytarabine, gemcitabine, cladribine and fludarabine) including nucleoside analogs of non-physiological L-chirality (e.g. 3TC, lamivudine, anti-HIV and anti-hepatitis B agent). Three-dimensional structures of dCK in complex with various pyrimidine<sup>1</sup> and purine<sup>2,3</sup> D- and L-nucleosides have been solved providing structural basis for activation of L- and D-nucleoside analogs.

**Storage:** -20 °C in a solution containing 50 mM Tris-HCl, pH 7.6, 1 mM β-mercaptoethanol, 50% glycerol.


**Unit Definition:** One unit of deoxycytidine kinase converts 1.0 μmole of deoxycytidine and ATP to dCMP and ADP per minute at pH 7.6 at 37°C, as measured by a coupled PK/LDH enzyme system.

**Specific Activity:** ≥ 0.025 unit/mg protein.

**Purity:** controlled by 12%AA SDS-PAGE.



The enzymatic activity of human recombinant dCK was confirmed by ion-pair HPLC analysis (Agilent 1100 series, Zorbax C18plus) as shown by formation of dCMP and ADP (red) from deoxycytidine and ATP (blue).



**Assay condition:** Enzymatic activity of dCK is measured by spectrophotometric assays in a coupled lactate dehydrogenase/pyruvate kinase system. Assays were carried out at 37°C, at 50mM Tris-HCl pH7,6; 50mM KCl, 10mM MgCl<sub>2</sub>, 5mM ATP, 0,1mM NADH, 1mM phosphoenolpyruvate, 1mM DTT, PK 10U/ml, LDH 15U/ml, 0,9μM dCK. Reaction was followed in an iEMS Reader MF (Labsystems) microtiter plate reader at 340nm. Nucleosides, nucleotides, LDH and PK were purchased from Sigma-Aldrich.

### Related products:

**NOVO CIB** has cloned and purified a panel of human recombinant nucleoside kinases and has developed a range of PRECICE® services to evaluate substrate properties of new nucleoside analogues for key cellular kinases.

- **dCK nucleoside phosphorylation assay**
- **Coupled dCK-CMK nucleoside phosphorylation assays**
- **Coupled Nucleoside Kinase – IMPDH II**
- **UMP-CMP kinase (CMK)**
- **Adenosine kinase (AK)**
- **Cytosolic 5' nucleotidase II (cN-II)**
- **CMK nucleotide monophosphate phosphorylation assay**
- **Adenosine kinase phosphorylation assay**
- **cN-II phosphorylation assay**

<sup>1</sup> Sabini E, Ort S, Monnerjahn C, Konrad M, Lavie A. **Structure of human dCK suggests strategies to improve anticancer and antiviral therapy.** (2003) *Nat Struct Biol.*10(7):513-9

<sup>2</sup> Sabini E, Hazra S, Ort S, Konrad M, Lavie A. **Structural basis for substrate promiscuity of dCK.** (2008) *J Mol Biol.* 378(3):607-21

<sup>3</sup> Sabini E, Hazra S, Konrad M, Burley SK, Lavie A. (2007) **Structural basis for activation of the therapeutic L-nucleoside analogs 3TC and troxacitabine by human deoxycytidine kinase.** *Nucleic Acids Res.* 35(1):186-92