



## PRECICE® Freshness Assay Kit

For rapid measurement of **IMP**, **inosine** and **hypoxanthine** in muscle tissues

Ref: 0700-003

Key words: K-value, ATP nucleotide catabolites

### I. Principle of PRECICE® Freshness Assay Kit

The "PRECICE® Freshness Assay" is a first microplate assay for assessment of quality of seafood freshness. The enzymes provided with kit allow specific conversion of three major ATP catabolites - IMP, inosine (Ino) and hypoxanthine (Hx) to NADH<sub>2</sub>. The quantification of nucleotides is done by measuring sample absorbance at 340nm. **PRECICE® K (IMP) Assay Kit** is based on the use of original recombinant enzymes of nucleotide metabolism that allow a simple and reliable quantification of **IMP**, **Ino** and **Hx** by measuring absorbance at 340nm (patented).

Reaction buffer + co-factors 1 and 2	Auto-absorbance of fish extract and co-factors	Blank
Enzyme mix 1	IMP → NADH	Absorbance <sub>340</sub> 1
Enzyme mix 2 + co-factor 3	IMP + Hx → NADH	Absorbance <sub>340</sub> 2
Enzyme Mix 3	IMP + Hx + Ino → NADH	Absorbance <sub>340</sub> 3

In each sample, **IMP** is totally converted to NADH<sub>2</sub> by specific dehydrogenases found in **Enzyme Mix 1**, while enzymes found in **Enzyme Mix 2** convert both **IMP** and **Hx**, and **Enzyme Mix 3** converts to NADH<sub>2</sub> all three major ATP catabolites - **IMP** + **Ino** + **Hx**. The NADH<sub>2</sub> formed is quantified by measuring specific absorbance at 340nm. The non-specific absorbance of sample value due to the presence of fish extract and co-factors is measured as a Blank and subtracted from experimental values. The absorbance in a long UV range (340nm) allows using common consumable and running the assay on a 96-well plate. For each sample, **IMP** (%), indicator of ultrafreshness, **Hx** (%) indicator of alteration and **Ino** (%) are calculated as follows:

$$\text{IMP (\%)} = \frac{\text{Absorbance}_{340}1 - \text{Blank}}{\text{Absorbance}_{340}3 - \text{Blank}} \times 100$$


$$\text{Hx (\%)} = \frac{\text{Absorbance}_{340}2 - \text{Absorbance}_{340}1}{\text{Absorbance}_{340}3 - \text{Blank}} \times 100$$

### II. Equipments required (not provided)

- 1) Boiling water bath
- 2) 50-ml tubes or bottles with screw caps resistant to heating (ex. polypropylene tubes from Corning ref. 430828)
- 3) 0.2µm non-sterile filters (ex. Sartorius)
- 4) 10mL non-sterile syringes (ex. Braun)
- 5) Plate agitator
- 6) Plate reader fitted with a 340nm filter (ex. Labsystems iEMS Reader MF (Thermo), Epoch (BioTec); PerkinElmer)



### III. OPTIONAL: Defining the linearity range of absorbance at 340nm of your instrument



In "PRECICE® K" kit, the quantification of ATP-degradation products totally relies on the NADH<sub>2</sub> absorbance. The linearity range of the instrument used for the quantification is of **critical importance** and should be checked before starting the experiments. NovoCIB provides upon request the standard solutions of NADH<sub>2</sub> that can be used for the calibration of spectrophotometer or plate reader.

We strongly recommend using plate reader because this instrument has a larger linearity range compared to a spectrophotometer and allows simultaneous reading of 96 samples without additional pipeting and liquid handling. Alternatively, the results can be read with a spectrophotometer by measuring absorbance at 340nm in disposable cuvettes (PMMA, 1.5ml ref. 759115, Brand). These cuvettes have to be filled at least with 400µL for absorbance reading, which is convenient since each 200µL assay must be diluted 2x or 4x before absorbance reading with a spectrophotometer.


### IV. Kit Content (one 96-well plate, 22 samples or 11 samples in duplicate)

- 1) "Extraction buffer concentrate", for preparing 1L of extraction solution.
- 2) "Reaction buffer" (provided in 15-ml tube)
- 3) "Enzyme mix 1", lyophilized in 15-ml glass vial, for IMP quantification
- 4) "Enzyme mix 2", lyophilized in 10-ml glass vial, for IMP + Hypoxanthine quantification
- 5) "Enzyme mix 3", lyophilized in 10-ml glass vial, for IMP + Hypoxanthine + Inosine quantification
- 6) 1 tube "Cofactor 1" (powder)
- 7) 1 tube "Cofactor 2" (powder)
- 8) 1 tube "Cofactor 3" (powder)
- 9) Transparent empty microplate (round-bottom 96-well plate Corning, Costar® ref. 3797) pre-filled with positive control (supplied dried)



## V. Storage

PRECICE® K-Freshness Assay Kit must be stored at -20°C until used. Enzyme mixes must be freshly prepared before performing the assays. Once prepared, reagent and enzyme solutions are not stable enough for storage.



## VI. Sample preparation

**VI.1. Before starting:** Thaw "Reaction buffer" tubes at room temperature (**1h in advance**). DO NOT HEAT! Since the rate of enzymatic reaction depends on the temperature, it is important to completely thaw the "Reaction buffer" and to equilibrate it at room temperature. Thaw "Extraction buffer".

**VI.2.** To prepare extraction buffer, quantitatively transfer the content of "Extraction buffer concentrate" vial into a bottle with 1L of deionized water.

**VI.3.** Weight 2-5g of muscles into a 50-ml polypropylene tube with screw caps resistant to heating (ex. polypropylene tubes from Corning ref. 430828), add 8-10\*volumes of extraction buffer per gram of muscle (see Table 1), tightly close the tube to avoid evaporation, and put them into a boiling water bath for 20min. Be sure that the tubes are put in the water deeply enough to cover all the muscle.



**Table 1.**

Raw material	Recommended dilution
Fish fillet	10 volumes of extraction buffer per gram
Shrimps	8 volumes per gram
Scallop	4 volumes per gram
Frog legs	1 volume per gram

Nucleotide concentration in muscle is ~5-10mM. The recommended dilutions allow to obtain ~300µM concentration of nucleotide in extract, that, once converted to NADH<sub>2</sub>, would correspond to less than 1.5 Absorbance Unit, within linearity range of common plate reader. This dilution also helps to avoid depletion of cofactors that would lead to incomplete enzymatic reactions and underestimated values.

**For additional information contact us directly at [contact@novocib.com](mailto:contact@novocib.com)**



**VI. 4. Filter the exudates:**


- Take off the plunger from syringe,
- Fix 0.2µm filter on the barrel
- Carefully transfer the exudate to the barrel
- Insert the plunger and filter the exudates into a clean tube.

**IMPORTANT: Do not filter hot or warm samples!**

**VII. Microplate filling with extracts (22 samples or 11 samples in duplicate per plate)**


**VII.1.** The microplate provided with the kit is pre-filled with a positive control (well H12 closed with a sticker). Peel off the sticker covering the well with the standard before use, add 100µL of “Extraction buffer” already dissolved in 1L water to all the wells of column 12 for standard solubilization and blanks.

100µL of Extraction buffer per well




**VII.2.** Dispense 100µL of diluted filtered extracts in a 96-well plate, as shown in the schemes.

**Scheme for 22 samples**



100µL of sample per well

## Scheme for 11 samples in duplicata



100µL of sample per well

## VIII. Preparation of reaction mixtures

**VIII. 1.** To solubilize lyophilized enzymes, add 200µL of deionized water to glass vial "Enzyme Mix 1" and 100µL of deionized water to 7-ml vial "Enzyme mix 2" Agitate gently until complete dissolution of the powder in the vials.

**VIII. 2.** Quantitatively transfer cofactors 1 and 2 to "Reaction buffer". To do so:


- pipet 1ml of "Reaction buffer" to each tubes and mix them by inverting or pipeting up and down until dissolved,
- transfer by pipeting the content of all two small tubes back into a 15-ml tube with "Reaction buffer",
- repeat to be sure that all reagent and enzymes of the small tubes are recovered.
- mix the complete "Reaction buffer" then containing cofactors 1 et 2 by gently inverting. Avoid bubbles.

**IMPORTANT: Do not add cofactor 3!**

**VIII. 3.** Pipet 9 mL of the so prepared "Reaction buffer" into the glass vial "Enzyme mix 1". Mix gently by inverting (do not vortex, avoid bubbles). "Enzyme mix 1" is ready.

**VIII. 4.** Transfer 6 mL of the "Enzyme mix 1" into the glass vial "Enzyme mix 2". Quantitatively transfer the content of the tube with cofactor 3 to "Enzyme mix 2". To do so, pipet 1ml of "Enzyme mix 2" to the tube with cofactor 3, mix by inverting. Avoid bubbles. "Enzyme mix 2" is ready.

**VIII. 5.** Transfer 3 mL of so prepared "Enzyme mix 2" into the vial "Enzyme mix 3" (5ml vial). Mix gently by inverting (do not vortex, avoid bubbles). "Enzyme mix 3" is ready.




You have:

- 3 ml of "Reaction mix" for Blank measurement
- 3 ml of "Enzyme mix 1" for IMP measurement
- 3 ml of "Enzyme mix 2" for IMP + hypoxanthine (Hx) measurement
- 3 ml of "Enzyme mix 3" for IMP + Inosine + hypoxanthine measurement



## IX. Microplate filling with reaction mixture

**IX. 1.** Dispense 100µL of "Reaction mix" (Blank) into the wells of line A and E; 100µL of "Enzyme mix 1" – of line B and F, 100µL of "Enzyme mix 2" - lines C and G, and 100µL of "Enzyme mix 3" – lines D and H as shown below.



**IX. 2.** Agitate for 2 min at 1000rpm. Incubate for 30min. Agitate again the plate before reading the absorbance at 340nm.

**IX. 3.** Program plate reader for single absorbance reading and read OD340 for whole plate. (Optional: You can follow Optical Density of the reaction at 340nm by programming plate reader for kinetics).

## X. Calculate nucleotide concentrations:

Copy-paste the values obtained with the single absorbance reading at 340nm (values after 30min incubation) in the Excel file, either in the 22 samples sheet or the 11 samples in duplicate sheet.

	Sample 1	Sample 3	Sample 5	Sample 7	Sample 9	Sample 11	Sample 13	Sample 15	Sample 17	Sample 19	Sample 21	Standards
Blank	0,251	0,246	0,240	0,239	0,242							0,263
Mix 1	1,044	1,030	0,491	0,944	0,777							0,272
Mix 2	1,111	1,081	0,570	1,020	0,828							0,3355
Mix 3	1,363	1,249	1,025	1,488	1,085							0,335
Blank	0,243	0,255	0,256	0,242	0,250							0,244
Mix 1	1,065	1,136	0,554	1,049	0,798							0,252
Mix 2	1,114	1,178	0,661	1,093	0,840							0,315
Mix 3	1,389	1,317	1,218	1,393	1,179							0,947
	Sample 2	Sample 4	Sample 6	Sample 8	Sample 10	Sample 12	Sample 14	Sample 16	Sample 18	Sample 20	Sample 22	

### References

- Saito, T, Arai, K.I. and Yajima T. (1959) Changes in purine nucleotides of red lateral muscle of rainbow trout. *Nature* 184: 1415
- Saito, T, Arai, K.I. and Tanaka, T. (1958) Changes in adenine nucleotides in squid muscle *Nature* 181: 127-1128
- Battle, N., Aristoy, M. C., & Toldrá, F. (2000). Early post-mortem detection of exudative pork meat based on nucleotide content. *Journal of Food Science*, 65(3), 413–416.
- Battle, N., Aristoy, M. C., & Toldrá, F. (2001). ATP metabolites during aging of exudative and nonexudative pork meats. *Journal of Food Science*, 66(1), 68–71.
- K. Saito a, A. Ahmed b, H. Takeda b, S. Kawahara b, M. Irie b, M. Muguruma Effects of a humidity-stabilizing sheet on the color and K-value of beef stored at cold temperatures. *Meat Science* 75 (2007) 265–272
- Luong, J.H.T., Male, K.B., Masson, C., & Nguyen, A.L. 1992. Hypoxanthine ratio determination in fish extract using capillary electrophoresis and immobilized enzymes. *J. Food Sci.*, 57, pp. 77 - 81.